



Fast Pulsed, High Current Laser Diode Driver

HIGHLIGHTS

- Fast Risetime 25ns
- High Current 200A (120A, 40A)
- Flexible Control Keypad or USB
- OEM version available

APPLICATIONS

- Lab
- R&D
- OEM

OPERATION PRINCIPLE

The LDD100 laser diode driver employs a digitally controlled analogue current source to provide for simple, flexible and accurate laser pulsing.

The unit may be USB-controlled for automated measurements. Alternatively, the user may set up all parameters and operate the unit via the front panel keypad and LCD display.

FIELDS OF APPLICATION

These laser drivers are particularly useful for laser diode testing procedures such as burn-in and LIV measurement. The LDD100 provides up to 200A in pulses ranging from 50ns to 3s duration. The rise time of the pulses is under 25ns with essentially no overshoot.

The flexibility in setup and control makes the unit very useful for general laboratory environments, where high power, pulsed diode lasers are used.



ABSOLUTE MAXIMUM RATINGS

Average Current	7A
Temperature Range	0 – 40 °C

ORDERING INFORMATION

Full order code: LDD 100 **c**

	Options	Description
Maximum current (c):	40	40A
	120	120A
	200	200A



For customized or OEM systems, please contact us.

SPECIFICATIONS


PARAMETER	CONDITIONS				UNITS
		MIN	TYP	MAX	
OUTPUT					
Current		400		200 ¹	mA A
Connector	Control unit to head Head to laser	LEMO Strip line			
Rise time			20	25	ns
Fall time			80	100	ns
Pulse duration		50		3	ns s
Resolution of pulse duration	50ns to 1.5ms 1.5ms to 3s	50		0.01	ns %
Repetition rate				100	kHz
Pulse separation		10		10	μs s
Resolution of pulse separation	10μs to 100ms 100ms to 1s 1s to 10s	10	100	1	μs μs ms
Compliance voltage				4	V
Duty cycle	without active cooling LDD100-40 LDD100-120 LDD100-200 with active cooling			0.012 0.004 0.003 0.8	
Accuracy		± 2			%
Reproducibility		± 1			%
Linearity			± 1		%
PC INTERFACE					
Type		USB			
Data transfer rate		1.5			Mbit/s
TRIGGER					
Type		Internal External: sequence start External: individual pulse control			
POWER SUPPLY					
Type		120V; 60Hz 240V; 50Hz			
DIMENSIONS					
	Control unit Current head	250 x 235 x 150 mm (W x L x H) 105 x 190 x 55 mm (W x L x H)			mm mm

¹ 40A and 120A versions also available

OPERATING MODES


Local via keypad and LCD display.

Remote via graphical user interface:




SEQUENCE MODES


Parallel staircase:



Serial staircase:



Pulse:



MENU STRUCTURE

In the following, the columns “Main menu” and the submenus show what appear on the LCD display in each line. The box outline indicates the full display in each case. The user can navigate these lines using the ↑ and ↓ keys. This results in rolling the menu up or down. Alternatively, the user can just press the number corresponding to the line and the menu is entered directly.

Main menu	Submenu 1	Submenu 2	Submenu 3	Submenu 4	
1/Sequence	1/Parallel	1/Current	1/Max	value + unit	
			2/Min	value + unit	
			3/Bias	value + unit	
		2/No. of steps	value		
			3/No. of cycles	Number	value
				Continuous mode	0=off / 1=on
	2/Serial	1/Current	1/Max	value + unit	
			2/Min	value + unit	
			3/Bias	value + unit	
		2/No. of steps	value		
			3/No. of cycles	Number	value
				Continuous mode	0=off / 1=on
3/Pulse	1/Current	value + unit			
	2/Bias	value + unit			
	3/No. of cycles	Number	value		
			Continuous mode	0=off / 1=on	
2/Trigger	1/Ext TRG: pulse sync				
	2/Ext. TRG: seq start				
	3/Int. TRG				
3/Pulse Timing	1/Width		value + unit		
	2/Period		value + unit		
	3/Separation		value + unit		
	4/Repetition rate		value + unit		
4/Extras	1/Current Limit		value + unit		
	2/Setup Profiles		1/Store 2/Recall		
			number (1-10)		

For example, on programme start (after the “intro”), the LCD displays:

1/Sequence
2/Trigger
3/Pulse Timing
4/Extras

The “Enter” key now enters a main menu and shows the available submenu. Each submenu is designed to have only 4 selections, so that no rolling is necessary or possible. For example, entering the main menu “Sequence” results in the following display:

1/Parallel
2/Serial
3/Pulse

From here, entering “1/Parallel” results in the following display:

1/Current
2/No. of steps
3/No. of cycles

From here, entering “1/Current” results in the following display:

1/Max
2/Min
3/Bias

From here, entering “1/Max” results in the following display:

<present value> <present units>

For example:

100000 mA

The user can enter a number and units here. The enter key then brings the display back to submenu 3.


In any given submenu, the “Cancel” button brings the display up one menu level.

The latest parameters are automatically saved into a “default” non-volatile memory when exiting a menu point or when starting a sequence. The user may begin a pulsing sequence from any main menu or submenu as long as that menu is not in edit mode. That is, as soon as the last parameter has been entered, the sequence may be started.


Additionally, the user may store or recall a setup profile into numerated non-volatile memories via the menu “Setup Profiles”. A total of 10 user defined profiles may be stored in this manner.

PERFORMANCE

OUTPUT PULSE SHAPE (40A PULSE)



SAFE OPERATING AREA



NOTICE

Artifex Engineering reserves the right to make changes to its products or to discontinue any product or service without notice and advises customers to obtain the latest version of relevant information to verify, before placing orders, that information being relied on is current and complete. All products are sold subject to the terms and conditions of sale supplied at the time of order acknowledgement, including those pertaining to warranty, patent infringement and limitation of liability.

Customers are responsible for their applications using Artifex Engineering components.

Artifex Engineering assumes no liability for applications assistance or customer product design. Artifex Engineering does not warrant or represent that any license, either express or implied, is granted under any patent right, copyright, or other intellectual property right of Artifex Engineering covering or relating to any combination, machine, or process in which such products or services might be or are used. Artifex Engineering's publication of information regarding any third party's products or services does not constitute Artifex Engineering's approval, warranty or endorsement thereof.

Attention: The LDD series of drivers is designed for use with lasers. Personnel who use this instrument must, therefore, be instructed in the safe use of lasers and laser beams.

Always wear the proper laser safety glasses designed for the laser in use!

Never allow the direct or reflected laser beam to impinge on the eyeball or to come into contact with the skin!

Copyright © 2010, Artifex Engineering



Dortmunder Str. 16-18
26723 Emden, Germany

Tel: +49-(0)4921-58908-0
Fax: +49-(0)4921-58908-29

info@afx-eng.com
<http://www.afx-eng.com>